

NH₃-Source

QUMA

Introduction:

The NH₃ source is an instrument to continuously produce pure NH₃ mixtures in humidified air with a flow rate of 0.5 - 2 l/min. NH₃ formation is based on the reaction of ammonium (NH₄⁺) with NaOH in diluted aqueous solution in a temperature controlled stripping coil similar to a HONO source developed by *Kleffmann et al.* showing fast time response of a few minutes (see Fig. 1). The NH₃ concentration is controlled by the known Henry's law constant and the instrument's parameter (NH₄⁺ and NaOH concentrations, gas and liquid flow rates, pH of the mixed solutions, temperature of the stripping coil).




Fig.1: Example of varying the NH₃ concentration by changing the NH₄⁺ concentration. Each data point reflects a 30 s period.

The time response of the NH₃-LOPAP instrument used here to measure NH₃ was 5 min.

The source produces pure NH₃/air mixtures humidified at the temperature of the stripping coil (dewpoint), which can be varied in the range 5-20°C. The NH₃ gas phase concentration (2-600 ppbv) is linearly correlated to the NH₄⁺ liquid phase concentration used and thus, the NH₃ concentrations can be easily adjusted (see Fig. 2). The precision of the source is typically ca. ±1 % (see Figure 3).




Fig. 2: NH₃ concentration of the NH₃ source as a function of the NH₄⁺ concentration used.

NH₃-Source

QUMA




Fig. 3: Precision of the NH₃ source. Noise may also be caused by the noise of the NH₃-LOPAP instrument used to measure NH₃ (1 % precision). Thus, the given precision of 1 % of the NH₃ source should be considered as an upper limit.

While the source produces stable and pure NH₃, the output of the source is depending on several variables (NH₄⁺ and NaOH concentrations, gas and liquid flow rates, pH of the mixed solutions, temperature of the stripping coil). The instrument is provided by a calculation sheet by which the output of the source can be absolutely calculated when these variables are measured. The estimated accuracy of this calculation is ca. 10 %, which is mainly limited by the uncertainties in the temperature dependent Henry's law constant of NH₃ in alkaline solutions. Thus, the source can be also used for approximate calibration of NH₃ instruments where simple calibration methods are not available.

References:

Kleffmann, J., T. Benter, P. Wiesen: Heterogeneous Reaction of Nitric Acid with Nitric Oxide on Glass Surfaces under Simulated Atmospheric Conditions, *J. Phys. Chem. A*, 2004, **108**, 5793-5799.